



LIGHTFORCE® DEEP TISSUE LASER THERAPY™

TREAT WITH CONFIDENCE

LightForce® Therapy Lasers are known around the world for innovative technology and evidence-based treatment solutions. With thousands of customers across the globe, Chattanooga® estimates over 1 million+ treatments are performed using LightForce® Therapy Lasers per month. LightForce® Therapy Lasers are unique amongst laser devices in having a singular focus in rehabilitation and pain management.

With over 7,000 articles published on Low Level Laser Therapy, photobiomodulation is well established in the scientific community for efficacy in treating musculoskeletal pain and increasing local blood circulation.

KEY FEATURES

- Empower IQ™ delivery system, smart handpiece, user friendly and responsive, durable long fiber (available with 25W XPi and 40W XLi).
- Exclusive LightForce® Treatment Heads with quick disconnect system.
- Exclusive massage ball technology allowing treatments on-contact.
- Influence® Technology treatment settings being adapted to specific patient data.
- High Power range from 15W to 40W reaching deep tissues and reducing treatment time.

INDICATIONS

Deep Tissue Laser Therapy™ is used to treat painful muscles and joints associated with acute and chronic soft tissue injuries. It is also indicated for helping relieve pain and stiffness associated with osteoarthritis via its ability to increase localized blood flow. These factors make therapeutic laser a sensible modality choice for post-activity recovery.



TECHNICAL INFORMATION

LASER TYPE: Class IV, Solid State

LASER WAVELENGTH: 980/810 nm

MINIMUM LASER POWER: **XLi/XPi** 1W **FXi** 0.5W

OPERATING MODES: CW or Pulsed

AIMING BEAM: 650 nm, < 4.5 mW

DIMENSIONS: **XLi/XPi** 17" x 10.5" x 11" (43 cm x 28 cm x 26 cm) **FXi** 15" x 11" x 8" (38 cm x 28 cm x 20 cm)

WEIGHT: **XLi/XPi** 19.6 lbs. (9 kg) **FXi** 7 lbs. (3,2 kg)

ORDERING INFORMATION

LIGHTFORCE® XLi



LIGHTFORCE® XPi



PART NUMBER	PART DESIGNATION DESCRIPTION
SYS-000656-000-00	LTS-4000-L-S8 (XLi) System - Int'l LightForce XLi (40W) Therapy Laser System

PART NUMBER	PART DESIGNATION DESCRIPTION
SYS-000837-000-00	LTS-2500-L-S8 (XPi) System - Int'l LightForce XPi (25W) Therapy Laser System

LIGHTFORCE® XLi STANDARD ACCESSORIES

KIT-000686-000-02	XLi SHP w/ Attachment Kit Case and Treatment Heads XLi Handpiece System Kit with 6 Treatment Heads (includes: 800um EmpowerIQ Delivery System, Attachment Kit Case, XL Rolling Pin, XL non-contact Cone, Large Contact Ball, Large Non-Contact Cone, Small Contact Ball, Flat Window)
PRD-000758-010-00-L	800um SHP v2 Assy Handpiece - 800um Empower IQ Delivery System
KIT-000685-000	XLi BDA Attachment Kit Case with 4 Treatment Heads Attachment Kit Case with 4 Treatment Heads (includes: Attachment Kit Case, Large Contact Ball, Large Non-Contact Cone, Small Contact Ball, Flat Window)
PRD-000229-010-00	BDA Attachment head Large Ball Alum - Rolling Treatment Head - Large contact Ball Attachment - Rolling
PRD-000230-010-L	BDA Attachment Head Large Cone Alum Treatment Head - Large non-contact Cone Attachment
PRD-000644-000	BDA Attachment Head XL Cone Treatment Head - XL Cone Attachment
PRD-000721-000	BDA Attachment Head XL Roller Treatment Head - XL Rolling Pin Attachment
PRD-000687-000	XL Band, LightForce, Black Black LightForce Band for XL Attachments
PRD-000683-000-01	LTS Carry Bag - LightForce Carry Bag for LightForce 25W / 40W

LIGHTFORCE® XPi STANDARD ACCESSORIES

KIT-000686-000-01	XPi SHP w/ Attachment Kit Case and Treatment Heads XPi Handpiece System Kit with 5 Treatment Heads (includes: 800um Empower IQ Delivery System, Attachment Kit Case, Large Contact Ball, Large Non-Contact Cone, Small Contact Ball, Small non-contact cone, Flat Window)
PRD-000758-010-00-L	800um SHP v2 Assy Handpiece - 800um Empower IQ Delivery System
KIT-000497-000-02	BDA Attachment Kit Case with 4 Treatment Heads (Aluminum Rolling) Attachment Kit Case with 4 Treatment Heads (includes: Attachment Kit Case, Large Contact Ball, Large Non-Contact Cone, Small Contact Ball, Small non-contact cone)
PRD-000229-010-00	BDA Attachment head Large Ball Alum - Rolling Treatment Head - Large contact Ball Attachment - Rolling
PRD-000230-010-L	BDA Attachment Head Large Cone Alum Treatment Head - Large non-contact Cone Attachment
PRD-000232-000-L	BDA Attachment head Small Cone Treatment Head - Small non-contact Cone Attachment
PRD-000683-000-01	LTS Carry Bag - LightForce Carry Bag for LightForce 25W / 40W

LIGHTFORCE® FXi



PART NUMBER	PART DESIGNATION DESCRIPTION
SYS-000829-000-00	LTC-1510-L-B6 (FXi w/ BHP) System - Int'l Light Force FXi (15W) Therapy Laser System

LIGHTFORCE® FXi STANDARD ACCESSORIES

PRD-000424-010-01	600um BHP Handpiece - 600um Empower Delivery System
KIT-000497-000-00-L	BDA Attachment Kit Case with 4 Treatment Heads Attachment Kit Case with 4 Treatment Heads (includes: Attachment Kit Case, LargeContact Ball, Large Non-Contact Cone, Small Contact Ball, Small non-contact cone)
PRD-000229-000-L	BDA Attachment head Large Ball Treatment Head - Large contact Ball Attachment
PRD-000230-000-L	BDA Attachment Head Large Cone Treatment Head - Large non-contact Cone Attachment
PRD-000232-000-L	BDA Attachment head Small Cone Treatment Head - Small non-contact Cone Attachment
ELX-000767-000-L	XP Power 150W 12V AC_DC Power Supply
PRT000297-L	Battery
PRD-000807-000-01	LTC Carry Bag - LightForce Carry Bag for LightForce 15W

STANDARD ACCESSORIES (ALL PRODUCTS)

PART NUMBER	PART DESIGNATION DESCRIPTION
PRD-000231-000-L	BDA Attachment head Small Ball Treatment Head - Small Contact Ball Attachment
PRD-000233-000	BDA Flat Window Treatment Head - Flat Window Attachment
PRT000240-L	Attachment Kit Case
PRT000269	LightForce Band, black Black LightForce Band for Large Attachments
PLS-000367-000-L	BDA, Accessory Dust Cap Attachment Dust protection Cap
PRT000374-L	Power Cord, "G" Style Power Cord, UK, "G" Type
PRT000103-L	Power Cord - EURO BLK 8FT (C Type) Power Cord, EU, "C" Type
CBL-000320-000-L	Power Cord, 2.9M BLK (I type China / Aus) Power Cord, Asia / Australia, "I" Type
PRT000011-L	Over Goggle Over Glasses Safety Goggles (OD 5+)
PRT000326	Regular Goggle, LightForce Green Regular Safety Goggles, LightForce Green (OD 5+)
PRT000473	P-Shield P-Shield Safety Goggles

OPTIONAL ACCESSORIES

PRD-000204-010-L	LTC Hard Rolling Case Hardside Rolling Case for 15W system
PRD-000215-010-L	LTS Hard Rolling Case Hardside Rolling Case for 25W / 40W system
PRD-000207-000-L	Laser Cart and Bin Cart / Bin

chattanooga®

lightForce®

CLINICALLY STUDIED TO
HELP IN THE TREATMENT OF:

Neck Pain ⁽¹⁾

Low Back Pain ^(2,3)

Sciatica ^(4,5)

TMJ ^(6,7)

Elbow and Joint Pain ⁽⁸⁾

Arthritis ⁽⁹⁾

1. Chow RT, Johnson MI, Lopes-Martins RA, Bjordal JM. Efficacy of low-level laser therapy in the management of neck pain: a systematic review and meta-analysis of randomised placebo or active-treatment controlled trials. *Lancet*. 2009 Dec 5;374(9705):1897-908. doi: 10.1016/S0140-6736(09)61522-1. Epub 2009 Nov 13. Erratum in: *Lancet*. 2010 Mar 13;375(9718):894. PMID: 19913903.
2. Alayat MS, Atya AM, Ali MM, Shosha TM. Long-term effect of high-intensity laser therapy in the treatment of patients with chronic low back pain: a randomized blinded placebo-controlled trial. *Lasers Med Sci*. 2014 May;29(3):1065-73. doi: 10.1007/s10103-013-1472-5. Epub 2013 Nov 2.
3. Yousefi-Nooraie R, Schonstein E, Heidari K, Rashidian A, Pennick V, Akbari-Kamrani M, Irani S, Shakiba B, Mortaz Hejri SA, Mortaz Hejri SO, Jonaidi A. Low level laser therapy for nonspecific low-back pain. *Cochrane Database Syst Rev*. 2008 Apr 16;(2):CD005107. doi: 10.1002/14651858.CD005107.pub4. PMID: 18425909.
4. Holanda VM, Chavantes MC, Silva DF, de Holanda CV, de Oliveira JO Jr, Wu X, Anders JJ. Photobiomodulation of the dorsal root ganglion for the treatment of low back pain: A pilot study. *Lasers Surg Med*. 2016 Sep;48(7):653-9. doi: 10.1002/lsm.22522. Epub 2016 May 2. PMID: 27135465.
5. Jovčić M, Konstantinović L, Lazović M, Jovčić V. Clinical and functional evaluation of patients with acute low back pain and radiculopathy treated with different energy doses of low level laser therapy. *Vojnosanit Pregl*. 2012 Aug;69(8):656-62. PMID: 22924260.
6. Abbasgholizadeh ZS, Evren B, Ozkan Y. Evaluation of the efficacy of different treatment modalities for painful temporomandibular disorders. *Int J Oral Maxillofac Surg*. 2020 May;49(5):628-635. doi: 10.1016/j.ijom.2019.08.010. Epub 2019 Sep 21. PMID: 31547949.
7. Khairnar S, Bhate K, S N SK, Kshirsagar K, Jagtap B, Kakodkar P. Comparative evaluation of low-level laser therapy and ultrasound heat therapy in reducing temporomandibular joint disorder pain. *J Dent Anesth Pain Med*. 2019 Oct;19(5):289-294. doi: 10.17245/jdapm.2019.19.5.289. Epub 2019 Oct 30. PMID: 31723669; PMCID: PMC6834715.
8. Roberts DB, Kruse RJ, Stoll SF. The effectiveness of therapeutic class IV (10 W) laser treatment for epicondylitis. *Lasers Surg Med*. 2013 Jul;45(5):311-7. doi: 10.1002/lsm.22140. Epub 2013 Jun 3. PMID: 23733499.
9. Dima R, Tieppo Francio V, Tower C, Davani S. Review of Literature on Low-level Laser Therapy Benefits for Nonpharmacological Pain Control in Chronic Pain and Osteoarthritis. *Altern Ther Health Med*. 2018 Sep;24(5):8-10. PMID: 28987080.



DJO Global
1a Guildford Business Park
Guildford, Surrey
GU2 8XG
UNITED KINGDOM

chattanoogarehab.com

© 2021 DJO - MKT-000-00011422-EN - REV A
Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.