

Chattanooga / Compex Professional Muscle Stim

Wireless Professional 4 Channel Full



The **Chattanooga Wireless Professional** is today's most advanced electrotherapy unit, offering enhanced ease-of-use and convenience while saving space and time to optimize patient treatment.

Muscle intelligence technology: Offers personalized stimulation, adapted to your patient's physiology.

By introducing the smart mi-Technology in each module and offering the best in market professional programs, this device reinvents functional rehabilitation as we know it:

A combination of e-stim with active exercise allowing you to improve therapy efficiency and patient satisfaction



'FUNCTIONAL REHABILITATION IN MOTION'

Wireless technology: Gives you and your patient total freedom of movement.

Remote control: Adapted to professional use with a silicon keypad which acts as a safeguard against water and gel penetration, and offers enhanced ergonomics.

Favourite programs: Direct access to your favourite programs; three clicks to start a session.



Docking station:

Offering module storage in the charging tablet and an extensible case to store electrodes and gels.

Travel facilities: For portable use, a modular system allows you to detach the charging tablet, and take it with you in a travelling case.

Remote control: Adapted to professional use with a silicon keypad which acts as a safe guard against water and gel penetration, and offers enhanced ergonomics.

Motor point pen: This useful tool together with the dedicated program helps ensure the optimal electrode placement.

Price: €1950+VAT

Dec Special Offer: €1650+VAT

Programs

Rehabilitation I: 5 programs

Disuse Atrophy. Reinforcement. Prevention of Disuse Atrophy. Muscle Lesion. Motor Point

Rehabilitation II: 12 programs

Hip prosthesis, Patellofemoral syndrome. ACL, Rotator Cuff, Lumbar Stabilization, Cardiac Rehabilitation. Ago/Antago Atrophy, Ago/Antago Reinforcement. Atrophy (Mod. Freq.), Reinforcement (Mod. Freq.). Atrophy for Haemophiliacs. Reinforcement for Haemophiliacs

Pain relief I: 6 programs

TENS (100Hz). Freq. Modulated TENS. Pulse Width Modulated TENS. Endorphinic. Burst TENS Alternated. Decontracting

Pain relief II: 12 programs

Knee Pain. Trapezius Pain. Shoulder Pain. Fracture Pain. Cervical Pain. Thoracic Back Pain. Low Back Pain. Lumbosciatica. Lumbago. Epicondylitis. Torticollis. Arthralgia

Vascular: 7 programs

Heavy Legs. Venous Insufficiency 1. Venous Insufficiency 2. Arterial Insufficiency 1. Arterial Insufficiency 2. Cramp Prevention. Capillarization

Neurologic rehabilitation: 4 programs

Hemiplegic Foot. Spasticity. Hemiplegic Shoulder. Slow Start Neuro Rehabilitation

Aesthetic: 6 programs

Toning. Firming. Shaping. Elasticity. Calorilysis. Adipostress

Conditioning I: 3 programs

Resistance. Strength. Active Recovery

Conditioning II: 12 programs

Potentiation. Endurance. Explosive Strength. Plyometry. Hypertrophy. Muscle Building. Low Back Reinforcement. Core Stabilization. Recovery Plus. Toning Massage. Relaxing Massage. Anti-Stress Massage

To Order: Tel: 062-77014

e-mail: sales@sportsphysio.ie

website: www.sportsphysio.ie

Chattanooga Theta

'FOR PROFESSIONALS WITH HIGH EXPECTATIONS'

THETA is a 4 channels stimulators developed to treat both standard (Pain, NMES) and specific protocols in rehabilitation, especially when triggered or co-contractions are requested. Neurology (Spasticity, Hemiplegic foot, Agonist-Antagonist) or specific pathologies like Post ACL Ligamentoplasty are protocols provided by the unit

MUSCLE INTELLIGENCE™ TECHNOLOGY

4 functions on all 4 channels: Mi- Scan, Mi-Tens, Mi-Range, Mi-Action (Active session with a voluntary contraction detected by the device to start the stimulation). 2 Mi-Cables included.

2+2

Choose one programme for channel 1 and 2, and then select another programme from the 2+2 programme list for channel 3 and 4. This allows treatment of two different body parts, or two patients simultaneously, or combines NMES and TENS as one treatment.

TRIGGER FUNCTION

Start the stimulation on demand by pressing a trigger button when needed (Hemiplegic foot, spasticity...)
Manage the session time and move directly to the next phase (recovery) when needed.

MOTOR POINT PEN

Very useful tool in combination with the dedicated programme helps ensure the optimal electrode placement.

PROGRAMMABILITY

Creates and personalizes programmes in order to customize the patient specific parameters.



Programs

Rehabilitation I: 5 programs

Disuse Atrophy. Reinforcement. Prevention of Disuse Atrophy. Muscle Lesion. Motor Point

Rehabilitation II: 12 programs

Hip prosthesis, Patellofemoral syndrome. ACL, Rotator Cuff, Lumbar Stabilization, Cardiac Rehabilitation. Ago/Antago Atrophy, Ago/Antago Reinforcement. Atrophy (Mod. Freq.), Reinforcement (Mod. Freq.). Atrophy for Haemophiliacs. Reinforcement for Haemophiliacs

Pain relief I: 6 programs

TENS (100Hz). Freq. Modulated TENS. Pulse Width Modulated TENS. Endorphinic. Burst TENS Alternated. Decontracting

Pain relief II: 12 programs

Knee Pain. Trapezius Pain. Shoulder Pain. Fracture Pain. Cervical Pain. Thoracic Back Pain. Low Back Pain. Lumbosciatica. Lumbago. Epicondylitis. Torticollis. Arthralgia

Vascular: 7 programs

Heavy Legs. Venous Insufficiency 1. Venous Insufficiency 2. Arterial Insufficiency 1. Arterial Insufficiency 2. Cramp Prevention. Capillarization

Neurologic rehabilitation: 4 programs

Hemiplegic Foot. Spasticity. Hemiplegic Shoulder. Slow Start
Neuro Rehabilitation

Conditioning I: 3 programs

Resistance. Strength. Active Recovery

Conditioning II: 12 programs

Potential. Endurance. Explosive Strength. Plyometry. Hypertrophy. Muscle Building. Low Back Reinforcement. Core Stabilization. Recovery Plus. Toning Massage. Relaxing Massage. Anti-Stress Massage

Price: €1390+VAT

Dec Special Offer: €1190+VAT

Chattanooga Rehab



Chattanooga Rehab

REHAB is a 4 channels stimulator developed for PTs and Rehab Clinics requesting a standard way of use of electro stimulation. It meets all of the requirements a professional demand: pain management, treatment of vascular diseases, and muscle stimulation with basic programs.

MUSCLE INTELLIGENCE™ TECHNOLOGY

4 functions on all 4 channels: Mi- Scan, Mi-Tens, Mi-Range, Mi-Action (Active session with a voluntary contraction detected by the device to start the stimulation). 2 Mi-Cables included.

2+2

Choose one programme for channel 1 and 2, and then select another programme from the 2+2 programme list for channel 3 and 4. This allows treatment of two different body parts, or two patients simultaneously, or combines NMES and TENS as one treatment.

TRIGGER FUNCTION

Start the stimulation on demand by pressing a trigger button when needed (Hemiplegic foot, spasticity...)

SKIP

Manage the session time and move directly to the next phase (recovery) when needed.

MOTOR POINT PEN

This very useful tool in combination with the dedicated programme helps to ensure the optimal electrode placement.

PROGRAMMABILITY

Creates and personalizes programmes in order to customize the patient specific parameters.

Price: €590+VAT

Dec Special Offer: €495+VAT

Programs

Rehabilitation I: 5 programs

Disuse Atrophy. Reinforcement. Prevention of Disuse Atrophy. Muscle Lesion. Motor Point

Pain relief I: 6 programs

TENS (100Hz). Freq. Modulated TENS. Pulse Width Modulated TENS. Endorphinic. Burst TENS Alternated. Decontracting

Vascular: 7 programs

Heavy Legs. Venous Insufficiency 1. Venous Insufficiency 2. Arterial Insufficiency 1. Arterial Insufficiency 2. Cramp Prevention. Capillarization

Conditioning I: 3 programs

Resistance. Strength. Active Recovery